

세로퀴크

Seroquick

HEPATITIS C ANTIBODY RAPID TEST DEVICE

INTENDED USE

Seroquick Hepatitis C Antibody Rapid Test Device is a rapid chromatographic immunoassay designed for the qualitative detection of antibodies to hepatitis C virus (HCV) in human serum, plasma, or whole blood samples. This test serves as a screening tool for identifying past or current hepatitis C virus infection, including in individuals who may be asymptomatic.

With its fast results and simple procedure, the test is suitable for use at home or in community settings, as well as in clinical and medical environments. It does not require specialized equipment, making it ideal for both well-equipped and resource-limited conditions. Early detection of HCV antibodies allows for timely medical evaluation and confirmatory testing, helping reduce the risk of liver complications and limiting transmission within high-risk and vulnerable populations.

PRINCIPLE

Seroquick Hepatitis C Antibody Rapid Test Device is based on a lateral flow immunochromatographic assay that employs a double-antigen sandwich detection method. The test strip contains recombinant hepatitis C antigens conjugated to colloidal gold particles, which are deposited onto the conjugate pad. When a sample of whole blood, serum, or plasma is applied to the sample well, it interacts with the gold-labeled antigens, and if anti-HCV antibodies are present, they form an antigen-antibody complex that migrates along the nitrocellulose membrane by capillary action.

As the complex moves upward, it encounters a second set of recombinant HCV antigens immobilized in the test region. These antigens capture the complex and produce a visible colored band, indicating the presence of antibodies to HCV. If the sample does not contain anti-HCV antibodies, no colored band appears in the test region. The presence of a control band confirms the validity of the test and proper sample flow, confirming that an adequate specimen was applied and fluid migration occurred as intended. The control band must be visible for the result to be considered valid. This internal procedural control ensures the reliability of the test even in decentralized or non-laboratory settings. The test delivers results within 15 minutes and requires no specialized instruments, making it suitable for home, clinical, and point-of-care use.

STORAGE

1. Store in a cool, dry place between 2°C and 30°C. Do not freeze.
2. Keep away from direct sunlight and excess humidity.
3. Do not open the sealed pouch until ready to use.
4. Use before the expiration date printed on the packaging.
5. Before use, allow all test components to gradually reach and stabilize at room temperature (15°C to 30°C) to ensure proper fluid migration.

6. For best results, test samples immediately after collection. Avoid using damaged or unsealed test devices.

PACKAGE CONTENTS

1. Instruction Manual
2. HCV Test Cassette
3. Buffer Tube with Dropper Tip
4. Automatic Sterile Lancet
5. Manual Sterile Lancet
6. Pipette
7. Alcohol Pad
8. Bandage

SPECIMEN COLLECTION

If the specimen is to be tested immediately, open the sealed foil pouch beforehand, and ensure that all test kit components are at room temperature and ready for use before specimen collection.

A. Whole Blood Specimen

1. Clean the fingertip with the alcohol pad and allow it to air dry.
2. Squeeze the fingertip firmly and pierce it using the sterile lancet.
3. Wipe away the first drop of blood using a fresh cotton pad.
4. Gently massage or squeeze the fingertip again to promote blood flow.
5. Use a disposable pipette to collect blood from the puncture site. Apply controlled pressure and avoid drawing the blood all the way into the bulb, as it may become trapped and difficult to dispense.

B. Plasma Specimen

1. Collect venous blood into an anticoagulant tube via venipuncture.
2. Mix the tube gently by inversion to prevent clotting.
3. Centrifuge the specimen at the appropriate speed and duration as per laboratory protocol to separate plasma.
4. Using a sterile pipette, carefully transfer the plasma into a properly labeled, sterile secondary tube, avoiding the buffy coat.

C. Serum Specimen

1. Collect venous blood into a plain or serum separator tube via venipuncture.
2. Allow the blood to clot at room temperature for 30 to 60 minutes.
3. Centrifuge the clotted blood as per protocol to separate the serum.
4. Carefully transfer the serum into a properly labeled, sterile secondary tube using a sterile pipette, avoiding any residual cellular material.


D. Specimen Handling and Storage

1. Test specimens as soon as possible after collection for optimal results.
2. If immediate testing is not feasible, store specimens at 2 to 8°C for up to 5 days.
3. For longer storage, freeze specimens at -20°C or lower.
4. Avoid multiple freeze-thaw cycles, as repeated freezing and thawing may compromise analyte stability.
5. Before testing, thaw frozen specimens slowly at room temperature and mix gently by inversion.
6. If specimens contain visible particulate matter, clarify by centrifugation prior to testing to ensure sample integrity.

TEST PROCEDURE

1. Open the foil pouch and remove the test device. If testing multiple specimens, clearly label the device with the corresponding patient or control identification. Place the device on a flat, clean, and dry surface.
2. Using the pipette, dispense 1 drop (50 µL) of whole blood, serum, or plasma into the sample well of the test cassette.
3. Dispense 2 drops (100 µL) of buffer solution to the same sample well.
4. Start the timer. Results may begin to appear as early as 30 seconds, but final interpretation should be made between 10 and 15 minutes.
5. Do not interpret results after 20 minutes, as results may become unreliable.

READING RESULT



Negative Result: If only the C band is visible and no band appears in the T region, this indicates that the specimen contains no detectable anti-HCV antibodies. The result should be interpreted as negative.

Positive Result: If a colored band appears in the T region along with the C band, this indicates the presence of anti-HCV antibodies in the specimen. The result is interpreted as positive.

Invalid Result: If the C band does not appear, the test is invalid regardless of a visible band in the T region. The test should be repeated using a new device.

Note: The intensity of the band in the T region may vary depending on the concentration of anti-HCV antibodies present in the specimen. A darker band indicates a higher antibody concentration. Even a faint band in the T region should be considered positive if the C band is also present.

PERFORMANCE

A. Sensitivity

Seroquick Hepatitis C Antibody Rapid Test Device was evaluated against a third-generation enzyme immunoassay (EIA), a well-established reference method for detecting anti-HCV antibodies. In this study, all 74 anti-HCV positive specimens confirmed by the reference method were correctly identified by Seroquick, resulting in a sensitivity of 100% within the evaluated sample set. Designed for both clinical and point-of-care use, Seroquick delivers rapid, visual results without requiring specialized laboratory equipment. These results demonstrate the test's high diagnostic accuracy and strong potential for effective hepatitis C screening across various settings.

B. Specificity

Seroquick Hepatitis C Antibody Rapid Test Device was evaluated for specificity using 221 specimens confirmed negative by a third-generation enzyme immunoassay (EIA). All 221 specimens were correctly identified as negative by Seroquick, yielding a calculated specificity of 100% within the evaluated sample set. This perfect specificity underscores the test's strong ability to accurately rule out non-infected individuals, making it highly dependable for rapid hepatitis C screening across a wide range of healthcare settings. Such diagnostic precision minimizes the risk of false-positive results and unnecessary follow-up testing, supporting confident decision-making in both clinical and point-of-care use.

C. Performance Summary

Parameter	Result	Notes
Total Samples	295	74 positive, 221 negative
Sensitivity	100% (74/74)	Verified by third-generation EIA
Specificity	100% (221/221)	All negatives correctly identified
Overall Accuracy	100% (295/295)	Correctly identified samples
True Positives	74	All results correctly identified
False Negatives	0	No missed infections
True Negatives	221	Verified by EIA
False Positives	0	No incorrect positives
Time to Result	~15 minutes	Typical visual result time

D. Cross-reactivity

A cross-reactivity study was conducted to assess the potential for analytical interference from substances commonly found in specimens from patients undergoing hepatitis C screening. The evaluation included a variety of prescription medications, recreational drugs, blood additives, endogenous compounds, and antibodies potentially present in co-infections or autoimmune conditions.

All substances were tested at concentrations exceeding normal therapeutic or physiological levels to simulate worst-case conditions. Under these test conditions, none of the evaluated substances caused false positive or false negative results in Seroquick Hepatitis C Antibody Rapid Test Device, supporting the test's specificity and robustness in diverse clinical populations.

Parameter	Conc.	Parameter	Conc.	Result
Interferon- α	10,000 IU/mL	Acetaminophen	200 μ g/mL	No interference
Methadone	1,000 μ g/mL	Ibuprofen	200 μ g/mL	No interference
Zidovudine (AZT)	400 μ g/mL	Lamivudine (3TC)	400 μ g/mL	No interference
Anti-HBc antibodies	500 IU/mL	Anti-HIV antibodies	500 IU/mL	No interference
Anti-dsDNA antibodies	500 IU/mL	Rheumatoid factor (RF)	500 IU/mL	No interference
Ethanol	1%	Methamphetamine	500 μ g/mL	No interference
Cocaine metabolite	300 μ g/mL	Ampicillin	200 μ g/mL	No interference
Albumin	5 mg/mL	Hemoglobin	2 mg/mL	No interference
Bilirubin (conjugated)	2 mg/mL	Glucose	2 mg/mL	No interference
Heparin	1%	EDTA	800 μ g/mL	No interference

LIMITATIONS

1. This product is intended for in vitro diagnostic use only and should not be used for any other purpose.
2. Environmental factors such as excessive humidity, extreme temperatures, or improper storage may adversely affect performance.
3. Interfering substances in the specimen, as well as technical or procedural errors, may lead to false results. These are factors that may fall outside the manufacturer's control.
4. While the test demonstrates high accuracy and performance comparable to laboratory-based methods, a small possibility of false positive or false negative results still exists.
5. In cases of inconclusive or unexpected results, further evaluation with confirmatory laboratory testing is recommended. As with all diagnostic tools, results should be interpreted alongside the patient's clinical presentation, medical history, and other relevant laboratory findings.



Seroquick Biotech Co., Ltd.
309-1 Seonin-ri, Sacheon-eup, Sacheon-si,
Gyeongsangnam-do, South Korea
seroquick.co.kr
Need help? Email us at en_help@seroquick.co.kr